CHARACTERISTIC OF STAINLESS STEELS —

TASK instruments are basically made from the martensitic stainless steel. It is the most suitable material known now for the dental instruments because the proper hardness, strength, and toughness can be realized by heat treatment. The stainless steel is a kind of metal alloy that forms on surface stable film effective to rustproofing but rust spreads quickly once the surface is damaged. Please care the instruments as following procedures to use for a long term.

USUAL CARE —

- Lubricate periodically on all instruments with hinges, locks and the other moving parts such as forceps, pliers, scissors and needle-holders.
- Sharpness of blades decreases by wear-out or damage. When the blades becomes dull, ask us for sharpening or reconditioning without delay.

CARE AFTER USE —**Washing : Flush visible stains and tissue**

Concise processing – washing by hand / ultrasonic cleaner
 Chemical processing – cleaning agent / ultrasonic cleaner with cleaning agent

- ① Check if the instrument is in good condition to use it again.
- ② Flush blood, fat or protein from the surface before it dries.
Washer disinfectors are not recommended to use.
- ③ Cleaning agent in which rust prevention is effective is recommended. Use it according to directions of the manufacturer.
- ④ Water temperature from 40°C to 60°C (104° F to 140° F) is the most effective in the ultrasonic cleaner. Place instruments in open position in a way that they do not touch each other on a tray in order to avoid insufficient cleaning.
- ⑤ Rinse the instruments thoroughly with water.
- ⑥ Inspect the instruments have been properly cleaned.

Disinfection : Making power of harmful microorganism lower to the human body

Concise processing – Boiling disinfection / Filtration
 Chemical processing – Disinfectant (chemical solution)

- ① Dry the instruments completely after sterilization. Oiling at the joint is recommended for the instruments with moving parts such as forceps, pliers, and needle-holders.
- ※ Glutaral is recommended to use. Do not use the following solution: Sodium hypochlorite, Benzalkonium chloride, Benzathonium chloride, Povidone-Iodine, Phenol, Formaldehyde, Chlorhexidine Gluconate, Functional Water (super-acid water, oxidation water), Household detergent

Sterilization : Removal of microorganism for asepsis

Concise processing – Autoclave
 Chemical processing – EOG (ethylene oxide gas)

- ① Recommended the sterilization by autoclave or by EOG.
 Dry Heat Sterilization, Chemiclave and Cold Sterization are not recommended.
- ② Use the autoclave below temperature of 132°C (269.6 ° F) in dry room.
- ③ Do not use the high pressure steam sterilization with the medicine made of alcoholic content.
- ④ Dry instruments completely against the risk of corrosion. The inner parts of hinges need more time to dry.
- ⑤ Store instruments in dry and dust-free places, and keep away from chemical products.

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This registered trade-mark is stamped on all our instruments made of titanium. Light, strong, rustless, titanium has the excellent properties which conform to biocompatibility and modern dental treatment.

※The above logo marking is printed on all titanium instruments.

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
The high grade stainless of “400 Series” has a proper hardness and withstands wear and corrosion. It is ideal for the dental instruments.

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The instruments are made of hard chrome steel and chrome plated for anticorrosive.

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
The pliers are inserted with a hardened steel tip material for a longer working life.


※Non inserted products are shown without the above circled illustration.

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
●2枚ジョイント TWO PIECES JOINT

●3枚ジョイント THREE PIECES JOINT (BOX JOINT)

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
The plier consists of two pieces which are not jointed by the main screw A directly but by intervention C. Since the removable main screw pressed by the screw B with a hexwrench from the opposite side, the joint is properly kept tightened and does not become loose.

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A light weight large diameter hollow handle minimizes the fatigue of the wrist.

※All instruments with the hollow handle are illustrated with the magnified view as shown above.